

FILLMED
LABORATOIRES



NCTF® BOOST 135 HA

PATIENT AFTERCARE

This brochure is designed for those patients who have been recommended NCTF® BOOST 135 HA by their healthcare practitioner.

Bio-Revitaliser

You have received this information because your healthcare provider has recommended NCTF® BOOST 135 HA for you.

To achieve optimal results from your treatment, it is advised to follow the post-treatment recommendations listed in this aftercare brochure.



TREATMENT RESULTS



Immediate Glow



Collagen Stimulation



Wrinkle Reduction



Improved Skin Tonicity

INFORMATION REGARDING APPLICATION

NCTF may be applied via microneedling or post-laser. Microneedling is a minimally invasive cosmetic procedure that creates micro-channels in the skin through very fine needles. Laser devices similarly open channels allowing NCTF to take effect on the deeper levels of your skin. Both these applications break the skin barrier, so you should take particular care of your skin. Please follow the aftercare instructions below, minimising the chance of skin infection and irritation.

24 HOURS POST-TREATMENT

- Avoid touching the treated area to minimise any risk of infection
- Wash your hands before you touch your face
- Do not wash your face or apply any skincare products for 4-6 hours
- Do not apply make-up for at least 24 hours
- Avoid exercise that causes sweating

2-3 DAYS POST-TREATMENT

- Avoid using any active skincare products (Exfoliants, Retinol and Vitamin C)
- Avoid swimming, sauna, and spa environments.

2-4 WEEKS POST-TREATMENT

- Avoid exposing the treated areas to sunlight and UV light
- Refrain from undergoing procedures like peels, laser treatments, cosmetic tattoos, and waxing.

IMPORTANT INFORMATION

Your healthcare provider would have mentioned that in rare cases, there may be side effects that require attention.

If you experience any prolonged tenderness, pain that needs to be managed with relief, unusual bruising, an allergic reaction, itching or rash, contact your treating healthcare practitioner immediately.

