

FILLMED COLLAGEN MASK TREATMENT- FOR VITALITY, TONE AND RENEWED ELASTICITY WITH HYDRATION

- 1 First cleanse with Cleansing oil
- 2 Second cleanse with micellar water
- 3 Use exfoliating cream on entire face, wet fingertips and with circular motions massage for three minutes, remove with warm compress. (Can use steam)
- 4 Apply Glycopeel mask with silicone brush and leave on face for ten minutes, remove with silicone brush and remove any residue with tepid compress/facial cleansing cloth
- 5 Wipe once with perfection solution
- 6 Fold the collagen sheet in half
- 7 Pre-cut the collagen sheet to leave a space for nose and lips
- 8 Fold in 3
- 9 Soak the collagen sheet using moist sponges or facial cleansing cloth
- 10 Once soaked, unfold collagen sheet
- 11 Apply the mask on the face
- 12 Sculpt with FILLMED silicone brush directly on the face using cryo-serum
- 13 Leave it to set on the face for 10 to 15 minutes and follow CRYOLED facial steps with CRYOLED or cold handpiece device for five minutes
- 14 Last 5 to 10 minutes do hand and arm massage with nourishing body balm
- 15 Remove collagen mask without rinsing
- 16 Apply HXR eye cream to eyes and lips
- 17 Apply B3 Recovery cream to face
- 18 Apply youth protect



BEAUTY IS AN ART, WE MADE IT A SCIENCE

FILLMED

A DIVISION OF THE CRYOMED GROUP IN AUSTRALIA AND NEW ZEALAND

1300 346 448 info@cryomed.com.au cryomed.com.au