

These are Guidelines

Protocol Superficial Wrinkles

For the face, forehead, neck, décolletage and back of hands

- 1-2 Regen BCT tubes (depending on area size)
- Nappage (Multiple transdermal superficial injections) with 27-32G mesotherapy needle, 35° to 45°angle
- Alternatively, with mesotherapy gun (e.g., U225 Needle Concept)
- Let the A-PRP excess dry on the skin
- Do not wash treated areas for the next 4 hours
- Frequency: 3 treatments at 1 month interval D0 / D30 / D60
- Follow up treatment: between M6 and M12 (or as required)

Protocol Deep wrinkles and Dark Circles

- 1-2 Regen BCT tubes + 1 Regen ATS tube (depending on area size)
- Prepare 1ml syringe with 0.1 ml of ATS serum + 0.9 ml A-PRP immediately prior to injection; repeat as required
- Injections in the reticular dermis with a 30G x 25 mm flexible needle, 1ml for each injection site (Using retrograde surgical technique)
- For the periorbital area, injections of small amounts (0.1 ml) around the orbit, avoiding the angular vein (high risk of thrombosis); massage the treated area
- Frequency: 3 treatments at 1 month interval D0 / D30 / D60
- Follow up treatment between M6 and M12 (or as required)

Protocol Mature Skin

- 1-2 Regen BCT tubes + 1 Regen ATS tube (depending on area size), combination of Nappage and injections (Using retrograde surgical technique)
- Nappage (with A-PRP only): multiple transdermal superficial injections with 27-32G mesotherapy needle, 35° to 45° angle, or with mesotherapy gun (e.g., U225, Needle Concept)
- Prepare 1 ml syringe with 0.1 ml of ATS serum + 0.9 ml A-PRP immediately prior to injection; repeat as required
- Injections in the reticular dermis with a 30G x 25 mm flexible needle, 1 ml for each injection site (Using retrograde surgical technique) repeat as required
- Let the A-PRP excess dry on the skin
- Do not wash treated areas for the next four hours
- Frequency: 3 treatments at 1 month interval D0 / D30 / D60
- Follow up treatment: between M6 and M12 (or as required)

Protocol Scars & Stretchmarks

- 1 Regen BCT tube
- Prepare 1 ml PRP in a 1 ml syringe
- Nappage around the scars (transdermal injections + multiple superficial injections in the papillary dermis)
- Injections with a 27-32G needle (35° to 45° angle) to stimulate the superficial layers of the dermis
- Deep injections point by point, underneath the scar with a needle 32G x 4 mm (angle 90°)
- For hypertrophic scars, debride the surface of scar before applying A-PRP (e.g., surgical resection, dermabrasion, laser, peeling etc.)
- Frequency: Will depend on severity of scars – Stretchmark that are still pink may regain the normal skin colour. If Silver/white will not re-pigment.

If required Always aspirate prior to injecting