Efficacy and Safety of High-Power Diode Laser at 1060 nm for Skin Rejuvenation Treatment

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¹Laser Institute of Photomedicine, Barcelona, Spain

²Cocoon Medical, Barcelona, Spain

KEY TAKEAWAY


- The 1060nm high-power diode laser used in the current case obtains results comparable to those of solid-state Nd:YAG while still keeping the advantages of diode laser-based systems.
- Photographs showed a progressive improvement of the wrinkles, nasogenian fold depth, skin roughness and appearance.

- The patient experienced no pain or discomfort and even described the procedure as pleasant. Aside from a transient erythema that remained for 20 minutes no side effects occurred.

- 3D fringe projection measurements of the wrinkle confirmed a 18.8% volume improvement, a 15.5% reduction of the maximum depth and a 22% average depth decrease of the wrinkle.

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
3D images and fringe projection measurements of the wrinkle. Before: A and C; after: B and D.

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- Nasogenian fold 3D measurements showed an 18.5% volume improvement, a 34.9% reduction of the maximum depth and a 12.6% average depth decrease, 3D測量顯示，虎紋體積減少了18.5%，最大深度減少了34.9%，平均深度減少了12.6%。

3D images and fringe projection measurements of the nasogenian fold. Before: A and C; after: B and D.

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