SMARTLUX

HIGH INTENSITY PHOTOTHERAPY DEVICE

MEDMIX

Total Medical Solution
Company

MEDMIX CO.,LTD

Total Medical Solution Company

INTRODUCTION

HISTORY

GLOBAL CHAIN


MEDIA

Established in 2009 by cooperating with a group of specialists and professionals in the medical field, MEDMIX is dedicated to provide the best solutions for our customers who want to have health and a beauty every day.

We ceaselessly invest in R&D and devote our times and passions to improve our performances for innovative and creative products yet sticking to the core values which are good quality, reasonable price, and easy & costless maintenance.

MEDMIX promises to do our best to contribute our society for safe, healthy, and beautiful life by providing brilliant and effective medical device solutions.

Establishment of MEDMIX

2009~2010

- Established Medmix Co.,Ltd
- Launched SMARTLUX
- Acquired KFDA, Medical CE, IS13485

2011~2012

Establishment of R&D Center

- Established R&D Center
- Venture Company Qualification
- Acquired FDA approval

2013~2014

Movement of Headquarters

- Move the head office to Songdo
- Won the Korea Export Award of Million Tower
- Launched SMARTLUX SLIM

2015~2016

Launching New Devices

- Acquired MHLW from Japanese Gov. for the first time in Korea
- Launched SMARTLUX MINI
- Launched SMART VIEWER
- Launched TRIPLE1470

MEDMIX CO.,LTD

Total Medical Solution Company

INTRODUCTION

HISTORY

GLOBAL CHAIN

MEDIA

North America

- United States of America

Europe

- Belgium
- Germany
- France
- Poland
- Czech
- Italy

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
“Essential Care after Plastic Surgery.
That’s SMARTLUX which is getting popular
for Post-Treatment.”

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Articles

“Significant Increase in the Population of
Plastic Surgery. Attention on Post-Care.”

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
GLOBAL CHAIN

MEDIA

"Most Powerful Phototherapy Device, SMARTLUX".

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Articles

"The Implementation of the Advanced Photo-therapeutic Apparatus, SMARTLUX"

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The Implementation of the Advanced Photo-therapeutic Apparatus, Showing the SMARTLUX

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Maximizing the Effect of the Irradiation Equipment with Powerful Energy Output

The principle of the irradiator treating equipment starts from the mitochondria within the skin absorbing light energy and activate it. The stimulated mitochondria are actively producing ATP, and the blood circulation increases and collagen and elastin is produced in abundance. Undergoing this process, injury or disease, the skin received treatment is gradually restored. LED Light irradiator having been used for treatment. Output energy values were lower compared to expected effects maintenance issues has been steadily emerging. The emergency of LLLT (Low Level Laser Therapy) which can improve the intensive treatment area without damaging the skin tissue, and can treat even skin wounds and the recovery after surgery, facilitate recovery after surgery and surgical procedure, and deal with even hair growth and obesity, attracted extraordinary attention of the medical community. SMARTLUX at the head of SMARTLUX premium model in 2010 gained confidence in the LED light irradiator market with powerful energy output.

The major media reported technology innovation of SMARTLUX in detail. By using high-intensity SLD chip not existing common LED chip in the limited space of 19x41cm and by placing 2450 chips, the sink rate of the energy made in England etc

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A large number of prominent doctors choose SMARTLUX because of its high intensity and effect. SMARTLUX is the medical device recommended by famous doctors who have used it.

Magazine

- ✓ Dr. Koh Ik-soo (Petit Novel Plastic Surgery)
- ✓ Dr. Park Seong-su (Bong Bong Plastic Surgery)
- ✓ Dr. Jeong Chan-u (Leaders Clinic)
- ✓ Dr. Kim Beom-jun (Chung-ang Univ. Clinic)

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
SMARTLUX has been introduced on the Korean news channel, MBC with an explanation that its light has efficacy in various skin diseases such as herpes zoster, acne and inflammation.

News
MBC

Study results have shown that LED light is effective for treatment of skin diseases.

SMARTLUX

High Intensity Phototherapy Device

INTRODUCTION

PRINCIPLE

CLINICAL INDICATION


SMARTLUX VERSION

SMARTLUX delivers healing power and energy of light composed of red, blue, yellow, green, and IR with high-intensity SLDs (Super-luminous diodes) developed by SMART Phototherapy technologies. From wound healing and pain relief to improvement of dermatological conditions, anti-aging treatment, a boost to skin's radiance , SMARTLUX effects exceptional results with its various wavelengths selection and wide variety of irradiation angle.

Scientific Background

- ❶ Light penetrates into subcutaneous tissue.
- ❷ Mitochondria absorb the photon light energy and they are energized.
- ❸ The stimulated mitochondria produce more ATP, which stimulates cells to be reproduced faster and function lively like younger cells.
- ❹ The super luminous light promotes cell wall exchange and stimulates the microcirculation of blood. By increasing cell reproduction and improving blood circulation, more collagen and elastin are produced, which leads to reduce wrinkles and decrease healing time. Skin will look younger, plumper, and healthier within 2-6 weeks.

How SMARTLUX works

SMARTLUX helps to proceed below skin recovery steps faster

(Process of Skin Recovery)

Nominal Wavelength(nm)	Cell Types/Action Level								
	Inflammation			➔	Proliferation		➔	Remodeling	All
	Mast	Neutro	Macro	Fibro	Endo	Fibro-Myo	Keratino		
590-595				+++					
SMARTLUX RED 630-670	++	+	++	+++	+++	+	+++		
790	++			++					
SMARTLUX IR 830	+++	+++	+++	+	++	+++	+++		
904	-		+/-	-					

*Mast : Mast cells / Neutro : Neutrophils / Macro : Macrophages / Fibro : Fibroblasts / Fibro-Myo : Fibroblast to myofibroblast transformation / Keratino : Keratinocytes

Source
Calderhead, R. Glen. *Lasers in Dermatology and Medicine*, Springer London. 2012, 254

SMARTLUX

High Intensity Phototherapy Device

INTRODUCTION

PRINCIPLE

CLINICAL INDICATION

SMARTLUX VERSION

LIGHT COLOR	WAVELENGTH	CLINICAL INDICATION	APPLICATION
BLUE 	420nm	Sterilization	Acne
YELLOW 	590nm	Blood Circulation	Bruise, Swelling, Post-Plastic Surgery Treatment, Vascular Lesion, Whitening
RED 	635nm	Cell Reproduction	Skin Recovering, Wound Healing, Hair and Scalp Treatment
IR 	830nm	Reduction of Inflammation	Pain Relief

SMARTLUX

High Intensity Phototherapy Device

INTRODUCTION

PRINCIPLE


CLINICAL INDICATION

SMARTLUX VERSION

SMARTLUX SLIM

BRAND TYPE	PREMIUM BRAND
COMPOSITION	HEAD PART / BODY PART
LAMP TYPE	INDIVIDUAL LAMP TYPE (A lamp has only one light)
HEAD OPTIONS	FX PLUS (RED+IR), FX (RED+IR), DUAL (YELLOW+IR)
NUMBER OF SLD	FX PLUS: 2,450EA (RED: 1,120EA + IR: 1,330EA) FX: 1,200EA (RED: 700EA + IR: 500EA) DUAL: 2,450EA (YELLOW: 1,120EA + IR: 1,330EA)
WAVELENGTH	FX PLUS: RED-635nm / IR-830nm FX: RED-635nm / IR-830nm DUAL: 590nm / IR-830nm
INTENSITY (MAX)	FX PLUS: 65mW-135mW/cm ² FX: 45mW-90mW/cm ² DUAL: 18mW-95mW/cm ²

FX PLUS(RED+IR)

FX (RED+IR)

DUAL(YELLOW+IR)

SMARTLUX

High Intensity Phototherapy Device

INTRODUCTION

PRINCIPLE





CLINICAL INDICATION

SMARTLUX VERSION

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BRAND TYPE	DIFFUSION BRAND
COMPOSITION	ALL IN ONE TYPE
LAMP TYPE	MULTIPLE LAMP TYPE (A lamp has three lights)
HEAD OPTIONS	HEAD1 (BLUE+RED+IR), HEAD2 (YELLOW+RED+IR), HEAD3 (BLUE+RED+YELLOW)
NUMBER OF SLD	HEAD1: 4,320EA (1,440EA per color) HEAD2: 4,320EA (1,440EA per color) HEAD3: 4,320EA (1,440EA per color)
WAVELENGTH	HEAD1: BLUE-420nm / RED-635nm / IR-830nm HEAD2: YEL-590nm / RED-635nm / IR-830nm HEAD3: BLUE-420nm / RED-635nm / YEL-590nm
INTENSITY (MAX)	HEAD1: B+R+IR = 28mW-90mW/cm ² HEAD2: Y+R+IR = 20mW-75mW/cm ² HEAD3: B+R+Y = 20mW-90mW/cm ²

<Rotative Head>

SMARTLUX - MAIN FEATURES

High Intensity Phototherapy Device

HIGH INTENSITY

RESEARCH PAPERS

DIVERGENCE ANGLE

VARIOUS OPTION

COOLING SYSTEM

RACK STRUCTURE



SMARTLUX has extremely high intensity with high technology.

SMARTLUX



vs

LOW QUALITY DEVICE



RED
(635nm)

IR
(830nm)


BLUE
(420nm)

See the strength at the point of specific wavelength which has clinical indication.

The strength at the point of other wavelength does not have any clinical indication. (Useless)

*Measured by Spectrometer

Stable wavelength

Unstable wavelength

SMARTLUX - MAIN FEATURES

High Intensity Phototherapy Device

HIGH INTENSITY

RESEARCH PAPERS

DIVERGENCE ANGLE

VARIOUS OPTION

COOLING SYSTEM


RACK STRUCTURE

Low quality devices are same as LED lighting or traffic light.

LOW QUALITY
DEVICE

=

Low quality devices can't give you any medical effect at all due to poor intensity of light.

Do you want to buy "Expensive LED Lighting" or "Medical Device"?

SMARTLUX - MAIN FEATURES

High Intensity Phototherapy Device

- HIGH INTENSITY
- DIVERGENCE ANGLE
- COOLING SYSTEM
- RESEARCH PAPERS
- VARIOUS OPTION
- RACK STRUCTURE

Multiple research papers have been published by SMARTLUX.
 → Reputable, Reliable, Validated Medical Device!

Efficacy of Light-Emitting Diode Photomodulation of Erythema After Fractional Carbon Dioxide Laser: A Pilot Study

IN Y. OTS, MD,* BREM J. KIM, MD,* MYUNG N. KIM, MD,* CHAN SENG E. KIM, MD†

BACKGROUND: The most common side effects of fractional carbon dioxide (CO₂) laser resurfacing are transient erythema and edema. Light-emitting diode (LED) devices increase fibroblast activity and hasten wound healing. The current study was designed to evaluate the efficacy of LED photomodulation in reducing post-fractional CO₂ laser resurfacing erythema.

OBJECTIVES: To evaluate the clinical efficacy of LED photomodulation in reducing post-fractional CO₂ laser resurfacing erythema.

MATERIALS AND METHODS: Randomly selected facial halves of 10 Korean (III-IV) were treated using a 635-nm wavelength LED array immediately after resurfacing. Each participant was subsequently treated with LED daily for 10 days. Photographs, subjective physician assessment, and chromometer erythema results, with clinical improvement assessed using a 5-point grading scale.

RESULTS: The post-laser erythema resolved faster on the experimental side. Improvements noted according to physician assessment and chromometer significantly improved between the two sides were first noted on day 4.

CONCLUSION: Treatment using a 635-nm wavelength LED array decreases post-fractional CO₂ laser treatment erythema.

The authors have indicated no significant interest with commercial supporters.

Fractional carbon dioxide (CO₂) laser resurfacing has become widely used in the treatment of photodamage and scarring because of its marked clinical efficacy and low post-treatment side-effect profile. Transient erythema, edema, and edema of the treated skin are the most common side effects associated with treatment, whereas severe or permanent complications are rare.¹ Nevertheless, despite the limited recovery period after fractional laser resurfacing, patients often report that the associated skin erythema and edema prevent them from continuing immediately with their regular activities of daily living.

Light-emitting diode (LED) photomodulation is a non-thermal technology using light, which has been shown to have a photomodulatory effect on the skin. LED therapy has been shown to have a photomodulatory effect on the skin, including erythema and edema. LED therapy has been shown to have a photomodulatory effect on the skin, including erythema and edema. LED therapy has been shown to have a photomodulatory effect on the skin, including erythema and edema.

*Department of Dermatology, College of Medicine, Chung-Ang University, Seoul, S. Korea
 †Emergency Medicine, College of Medicine, Chung-Ang University, Seoul, South Korea

Correspondence

Efficacy of combination light-emitting diode (635 and 830 nm) therapy in treating local injection-site reactions after filler

doi: 10.1111/ced.12480

Over the past decade, the use of injectable cosmetic fillers has steadily increased, with hyaluronic acid (HA) preparations being the most frequently utilized. Complication rates with injections of HA fillers have been reported to be up to 5%. The most common reaction to HA is a local injection-site reaction, which consists of swelling, erythema and induration at the injection site.¹ Although these reactions are predominantly mild or moderate in intensity, and do not result in severe sequelae, non-treatment methods promoting rapid healing with high patient compliance are highly desirable because of safety concerns regarding the use of topical steroid agents as anti-inflammatory agents.² Consequently, researchers have sought other treatment options. Light and laser therapies have attracted attention as promising non-treatment methods for this condition.³

Light-emitting diode (LED) photomodulation is a non-thermal technology used to modulate cellular activity with light. Photomodulation refers to a photomodulatory effect using non-thermal cellular stimulation at specific pulse sequences and durations.⁴ Previous animal and human studies have demonstrated that LED treatment induces reduction in wound size and promotes wound healing. LED therapy using a variety of red, blue, and yellow wavelengths has been reported to accelerate cutaneous wound healing after various injuries, including inflammatory situations.⁵ To date, no studies have been published on the specific modulation of the efficacy of LED for the treatment of local injection-site reactions after filler. In the present study, we evaluated a combination LED (635 and 830 nm) therapy for treatment.

A 32-year-old woman received an injection of cross-linked HA filler (Kotylease Perlane®, QMED, Uppsala, Sweden) into the glabellar area in an attempt to plump up the area. The HA concentration in the filler was 20 mg/mL, and the total volume used was 1 mL. Two days after the injection, the patient returned with an inflammatory reaction in the injection area with an

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A Case of Small Plaque Psoriasis Treated with Combination of 578-nm Copper Bromide Laser with Light-emitting Diode (LED)

Eun Jung Ko, M.D., In Young Oh, M.D., In Pyeong S. M.D., So Young Kim, M.S., Tae Rin Kwon, M.S., Hye Mi Kim, M.S., Kwang Ho Yoo, M.D., Beom Joon Kim, M.D., Ph.D., Myeung Nam Kim, M.D., Ph.D.

Department of Dermatology, Chung-Ang University College of Medicine, Seoul, Korea

Psoriasis is characterized by chronic recurrent erythematous skin plaques that show enlargement and increased vascularity of cutaneous microvessels without formation of vessel sprouts, that is, inflammatory angiogenesis. The 578-nm copper bromide laser is effective modality for the vasculature. Light-emitting diode (LED) induce diminution in vessel size and protection from skin inflammation and ulceration report a case of small plaque psoriasis in 28-year-old successfully treated with combination of 578-nm copper bromide laser with 830-nm LED phototherapy. Clinical improvement was noted after 6 times treatment per week.

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REDUCING ERYTHEMA AFTER FRACTIONAL LASER

TREATMENT FOR LOCAL INJECTION-SITE REACTIONS AFTER FILLER

TREATMENT FOR SMALL PLAQUE PSORIASIS

RAPID REDUCTION OF POST-LASER RESURFACING PAIN, EDEMA, ERYTHEMA

SCI CLASS

SMARTLUX - MAIN FEATURES

High Intensity Phototherapy Device

HIGH INTENSITY

RESEARCH PAPERS

DIVERGENCE ANGLE

VARIOUS OPTION


COOLING SYSTEM

RACK STRUCTURE


SMARTLUX has a narrow divergence angle for intensified treatment result.

Divergence Angle

Divergence Angle

Narrow and precise positioning of the SLDs in the panel array results in multiple intersections of the light beams, generating the highest photon intensity zone on the target surface of the skin, providing intensified effects.

SMARTLUX - MAIN FEATURES

High Intensity Phototherapy Device

HIGH INTENSITY

RESEARCH PAPERS

DIVERGENCE ANGLE

VARIOUS OPTION

COOLING SYSTEM

RACK STRUCTURE

SMARTLUX provides a wide selection to meet various demands.

SMARTLUX SLIM **FX PLUS (R+IR)**

SMARTLUX SLIM **FX (R+IR)**

SMARTLUX SLIM **DUAL (Y+IR)**

SMARTLUX MINI **HEAD1 (B+R+IR)**

SMARTLUX MINI **HEAD2 (Y+R+IR)**

SMARTLUX MINI **HEAD3 (B+R+Y)**

SMARTLUX - MAIN FEATURES

High Intensity Phototherapy Device

HIGH INTENSITY

RESEARCH PAPERS

DIVERGENCE ANGLE

VARIOUS OPTION

COOLING SYSTEM

RACK STRUCTURE

SMARTLUX Cooling System lengthens device life time.

Temperature Controlling System

FAN


Cooling fans reduce the heats.

RADIATOR

Radiators lower the temperature.

SENSORS

Sensors maintain the temperature under 60°C

IN CASE OF USING SMARTLUX IN 50°C

IN CASE OF USING SMARTLUX IN 75°C

IN CASE OF USING SMARTLUX IN 125°C

SMARTLUX can keep its intensity at least for 10,000hours due to the cooling system which maintains the temperature under 60°C.

SMARTLUX - MAIN FEATURES

High Intensity Phototherapy Device

HIGH INTENSITY

RESEARCH PAPERS

DIVERGENCE ANGLE


VARIOUS OPTION

COOLING SYSTEM

RACK STRUCTURE


SMARTLUX has a rack structure for easy maintenance.

SMARTLUX consists of multiple racks.

Therefore, it's possible to repair the machine by just replacing the defective rack or component. It's very easy to manage the clients regardless of the distance between the distributor and clients.

No need to visit your customer's office or clinic to fix the machine.

If there is some problem, you Just send the defective rack or component to save time and money!

SMARTLUX – COMPARISON

High Intensity Phototherapy Device

KOREA - HEALITE II

CHINA - KN-7000A

FRANCE - TRIWINGS

SMARTLUX

Company	Medmix Co.,Ltd (Korea)
Number of SLDs	<ul style="list-style-type: none"> •SMARTLUX SLIM(FX PLUS): 2,450 EA •SMARTLUX MINI: 4,320 EA
Simultaneous Irradiation	<p>Possible You can use two or three colors at once.</p>
Wavelength Options & Intensity(Max)	<p>[SLIM]</p> <ul style="list-style-type: none"> •FX PLUS: R(635nm)+IR(830nm) → 65mW-135mW/cm² •FX: R(635nm)+IR(830nm) → 45mW-90mW/cm² •DUAL: Y(590nm)+IR(830nm) → 18mW-95mW/cm² <p>[MINI]</p> <ul style="list-style-type: none"> •HEAD1: B(420nm)+R(635nm)+IR(830nm) → 28mW-90mW/cm² •HEAD2: Y(590nm)+R(635nm)+IR(830nm) → 20mW-75mW/cm² •HEAD3: B(420nm)+R(635nm)+Y(590nm) → 20mW-90mW/cm²

HEALITE II

Company	Lutronic Co.,Ltd (Korea)
Number of LEDs	<ul style="list-style-type: none"> •1,800 EA
Simultaneous Irradiation	<p>Impossible You can't use two or three colors at once.</p>
Wavelength Options & Intensity(Max)	<p>[HEALITE II]</p> <ul style="list-style-type: none"> •IR(830nm) w/Y(590nm) → 40mW-100mW/cm² •R(633nm) → 30mW-65mW/cm² •B(415nm) → 10mW-30mW/cm² <p>[HEALITE II COMBO]</p> <ul style="list-style-type: none"> •R(633nm)+IR(830nm) w/Y(590nm) → 20mW-80mW/cm² •B(415nm)+R(635nm) → 10mW-50mW/cm² •B(415nm)+IR(830nm) w/Y(590nm) → 10mW-80mW/cm²

SMARTLUX – COMPARISON

High Intensity Phototherapy Device

KOREA - HEALITE II


CHINA - KN-7000A

FRANCE - TRIWINGS


SMARTLUX SLIM DUAL

YELLOW(590nm) + IR(830nm)

HEALITE II

YELLOW(590nm) + IR(830nm)

SMARTLUX SLIM DUAL

Company	Medmix Co.,Ltd (Korea)
Number of SLDs	•SMARTLUX SLIM(DUAL): 2,450 EA
Intensity(Max)	•Testing Result: 101mW/cm² (Official Data: 95mW/cm ²)
Yellow(590nm)	There is a distinct signal in 590nm.

HEALITE II 830(w/590)

Company	Lutronic Co.,Ltd (Korea)
Number of LEDs	•HEALITE II: 1,800 EA
Intensity(Max)	•Testing Result: 44mW/cm² (Official Data: 100mW/cm ²)
Yellow(590nm)	There is no signal in 590nm at all.

◇ Intensity is based on the actual survey.

SMARTLUX – COMPARISON

High Intensity Phototherapy Device

KOREA - HEALITE II

CHINA - KN-7000A

FRANCE - TRIWINGS

SMARTLUX

KN-7000A

Company	Medmix Co.,Ltd (Korea)
Number of SLDs	<ul style="list-style-type: none"> •SMARTLUX SLIM(FX PLUS): 2,450 EA •SMARTLUX MINI: 4,320 EA
Intensity(Max)	<ul style="list-style-type: none"> •SMARTLUX SLIM(FX PLUS): 135mW/cm² •SMARTLUX MINI(HEAD1): 90mW/cm²
Maximum Simultaneous Irradiation	3 Colors(Wavelengths) [MINI] <ul style="list-style-type: none"> •HEAD1: B(420nm)+R(635nm)+IR(830nm) •HEAD2: Y(590nm)+R(635nm)+IR(830nm) •HEAD3: B(420nm)+R(635nm)+Y(590nm)

Company	Kernel Co.,Ltd (China)
Number of LEDs	<ul style="list-style-type: none"> •KN-7000A: 1,820 EA
Intensity(Max)	<ul style="list-style-type: none"> •KN-7000A: 120mW/cm²
Maximum Simultaneous Irradiation	2 Colors(Wavelengths) <ul style="list-style-type: none"> •B(415nm)+R(635nm) •Y(590nm)+R(635nm) •R(635nm)+IR(830nm)

SMARTLUX – COMPARISON

High Intensity Phototherapy Device

KOREA - HEALITE II

CHINA - KN-7000A

FRANCE - TRIWINGS

SMARTLUX

Company	Medmix Co.,Ltd (Korea)
Intensity(Max)	<ul style="list-style-type: none"> •SMARTLUX SLIM(FX PLUS): 135mW/cm² •SMARTLUX MINI(HEAD1): 90mW/cm²
Energy Range	<ul style="list-style-type: none"> •SMARTLUX SLIM(FX): 1~324 J/cm² •SMARTLUX MINI(HEAD1): 1~184 J/cm²

TRIWINGS CLASSIC

Company	Biophoton Co.,Ltd (France)
Intensity(Max)	<ul style="list-style-type: none"> •TRIWINGS CLASSIC: 105mw/cm²
Energy Range	<ul style="list-style-type: none"> •TRIWINGS CLASSIC: 1~99 J/cm²

SMARTLUX - CLINICAL RESEARCH

High Intensity Phototherapy Device

ACNE

POST-LASER TREATMENT

WOUND HEALING

BRUISE & EDEMA

ECZEMA

FRECKLE

BEFORE

AFTER

SMARTLUX TREATMENT

SLD LIGHT	BLUE
LEVEL	4~5
TIME	10 MIN
FREQUENCY	TWICE A WEEK
PERIOD	2 WEEKS

SMARTLUX - CLINICAL RESEARCH

High Intensity Phototherapy Device

ACNE

POST-LASER TREATMENT

WOUND HEALING

BRUISE & EDEMA

ECZEMA


FRECKLE


BEFORE

1ST DAY

3RD DAY

SMARTLUX TREATMENT

SLD LIGHT	FX (RED+IR)
LEVEL	5
TIME	20 MIN
FREQUENCY	EVERYDAY
PERIOD	3 DAYS

SMARTLUX - CLINICAL RESEARCH

High Intensity Phototherapy Device

ACNE

POST-LASER TREATMENT

WOUND HEALING

BRUISE & EDEMA

ECZEMA

FRECKLE

BEFORE

AFTER

SMARTLUX TREATMENT

SLD LIGHT	FX (RED+IR)
LEVEL	5
TIME	20 MIN
FREQUENCY	2 TIMES
PERIOD	3 DAYS

SMARTLUX - CLINICAL RESEARCH

High Intensity Phototherapy Device

ACNE

POST-LASER TREATMENT

WOUND HEALING

BRUISE & EDEMA

ECZEMA

FRECKLE

BEFORE

AFTER

SMARTLUX TREATMENT

SLD LIGHT	FX (RED+IR)
LEVEL	5
TIME	20 MIN
FREQUENCY	1 TIME
PERIOD	1 DAY

SMARTLUX - CLINICAL RESEARCH

High Intensity Phototherapy Device

ACNE

POST-LASER TREATMENT

WOUND HEALING

BRUISE & EDEMA

ECZEMA

FRECKLE

BEFORE

1ST DAY

3RD DAY

SMARTLUX TREATMENT

SLD LIGHT	RED
LEVEL	5
TIME	20 MIN
FREQUENCY	EVERYDAY
PERIOD	2 DAYS

SMARTLUX - CLINICAL RESEARCH

High Intensity Phototherapy Device

ACNE

POST-LASER TREATMENT

WOUND HEALING

BRUISE & EDEMA

ECZEMA

FRECKLE

BEFORE

AFTER

SMARTLUX TREATMENT

SLD LIGHT	YELLOW
LEVEL	5
TIME	20 MIN
FREQUENCY	TWICE A WEEK
PERIOD	3 MONTHS

Thank you

sales@medmix.co.kr

MEDMIX

Total Medical Solution
Company