



CLATUU 













- Cooling Cup Features & Clinical Outcomes

- Table of Contents -

- 1. Cooling Cup Features & Target Areas**
 - 2. Clinical Outcomes**
 - 3. Clinical Background & Combinatory Treatments**
 - 4. Cooling Cup Features: a Comparison**
 - 5. Various Patient Postures During Treatment**
-

1. Cooling Cup Features & Target Areas

• Cooling Cup Features & Target Areas

 <p>• Wing</p>	 <p>• Wing Plus</p>	 <p>• Flat</p>	 <p>• Flat mini</p>	 <p>• Wing mini</p>	 <p>• Arc-T</p>	 <p>• Wing max</p>
<p>-For Large Fat pockets -Extending Applicator -For Broad & Curving Contours</p>	<p>-For Small Fat Pockets -Blunt-end Applicator -For Flat Contours</p>	<p>-For Large Fat Pockets -Flat-end Applicator -For Flat Contours</p>	<p>-For Small Fat Pockets -Flat-end Applicator -For Flat Contours</p>	<p>-For Small Fat Pockets -Applicable for Local Regions -For Broad & Curving Contours</p>	<p>-For Thigh Region -Wide area at once</p>	<p>-For Large Fat Pockets -Wide abdomen area at once</p>
<p>Ex) Upper Abdomen, Lower Abdomen, other large areas</p>	<p>Ex) Bra Line, Love Handles, Back</p>	<p>Ex) Abdomen, Inner Thighs</p>	<p>Ex) Double Chin (Side) Knees (Side) Other Local Regions</p>	<p>Ex) Axilla (Front & Back) Above the Knees</p>	<p>Ex) Around the thighs, flat and wide areas</p>	<p>Ex) Lower Abdomen</p>
 	 	 	 	 	 	

2. Clinical Outcomes

2. Clinical Outcomes - Wing



• Wing

- For Large Fat pockets
- Extending Applicator
- For Broad & Curving Contours

Ex) Upper Abdomen,
Lower Abdomen,
other large areas



BASELINE



4 Weeks Post Procedure

- Wide abdominal circumference with large fat pockets.
- Treats various target areas that are symmetrical.
- Application: Abdomen, flank, other fat accumulated areas

2. Clinical Outcomes - Wing Plus



• Wing Plus

- For Small Fat Pockets
- Blunt-end Applicator
- For Flat Contours

Ex) Bra Line, Love Handles, Below the Buttocks, Back



BASELINE



6 Weeks Post Procedure

- **Treats smaller fat pockets on bony and curving target area(s).**
- Designed for treatments particularly on the back with varying degrees of curvature.
- **Applications: Bra Line (Bra-fat), Back, Love Handles, Other.**

2. Clinical Outcomes - Flat



• Flat

- For Large Fat Pockets
- Flat-end Applicator
- For Flat Contours

Ex) Abdomen, Inner Thighs



BASELINE



4 Weeks Post Procedure

- Treatments applied on target areas with a flat surface.
- For broad areas with a concentrated fat pockets.
- **Applications: Abdomen and Inner Thighs.**

2. Clinical Outcomes - Wing mini



• Wing mini

- For Small Fat Pockets
- Applicable for Local Regions
- For Broad & Curving Contours
- Ex) Double Chin, Axilla (Front & Back) Above the Knees



BASELINE



4 Weeks Post Procedure

- For small, curving areas.
- More effective procedures and outcomes with the use of a cushion.
- Applications: Double Chin, Above the Knees, Axilla, Other.

2. Clinical Outcomes - Flat mini



- Flat mini

- For Small Fat Pockets
- Flat-end Applicator
- For Flat Contours

Ex) Double Chin (Lateral) Knees (inner) Other Local Regions



BASELINE



4 Weeks Post Procedure

- For small, flat areas.
- More effective procedures and outcomes with the use of a cushion.
- Applications: Double Chin (Lateral), Knees (Side), and other local regions.

2. Clinical Outcomes - Arc-T



• Arc-T

-For thighs
-Covering wider areas

Ex) Around the thighs, flat and wide areas



BASELINE



4 Weeks Post Procedure

- Designed for outer and inner thighs and below the buttocks.
- For effective suction, adjust the cup to where the target area is broadest.
- **One-time treatment to effectively remove excessive fat on thighs**

2. Clinical Outcomes – WING Max



- WING MAX

- Large area
- Bigger area at once
- Lower belly



BASELINE



4 Weeks Post Procedure


- Designed for covering bigger and wider areas of abdomen
- For effective suction, adjust the cup to where the target area is broadest and fat is mostly accumulated.
- **One-time treatment on broad areas of focus.**

3. Clinical Background & Combinatory


Treatments

- Double Chin
 - Abdomen
 - Back
 - Thighs
-

Clinical Background: Treatment on the Double Chin








- The subject may be suitable for treatment if excess fat accumulates in the submental triangle.



- Avoiding areas like the cheek bone, it is recommended to treat the sagging platysma muscle, loose skin, and neck with concentration of fat.
-

Clinical Background: Treatment on the Double Chin

- Eligible Subjects for Treatment (Scale of 2-4)

Scale	0	1	2	3	4
Submental Convexity	Absent	Mild	Moderate	Severe	Extreme
Description	No localized submental fat evident	Minimal localized submental fat	Prominent localized submental fat	Marked localized submental fat	Extreme submental convexity
Representative Photographs					

Combinatory Treatments – Double Chin



Double Chin (Side)



Double Chin (Side)




BASELINE



4 Weeks Post Procedure

Clinical Background: Treatment on the Abdomen


- **Devise a treatment plan according to body shape & profile**
 - Consider the direction/position of fat and muscle shape for the most appropriate treatment possible.



- **Subcutaneous vs. Visceral Fat**

- Excluding fat that is stored in the abdominal cavity (visceral fat), fat decomposition from exposure to cooling energy applies only on the subcutaneous fat region.

Combinatory Treatments – Abdomen



Abdomen (Curving Portion) + Abdomen (Flat Portion)



BASELINE



4 Weeks Post Procedure

Clinical Background – Treatment on the Back

- **Fat stretching from the spine to the body's outer edges is relatively short.**

- Fat in the back apportioned varies according to the location of fat accumulation.




- **Understanding the tissue profile and anatomy on the back for suitable cooling cup selection**

- Considering placement and direction of fat + obstruction (bones and ligaments)

- **Select the cup according to body shape and profile**

- In case the cooling cup is too lengthy or larger than the intended target area, breast tissue may also be suctioned into the cup. To prevent this, use the Wing Plus cooling cup instead of using Wing cooling cup.

Combinatory Treatments - Back



Back (Long, Curving Portion)

Back (Narrow, Gradual Portion)




BASELINE



4 Weeks Post Prodedure

Clinical Background – Treatment on the Thighs

- Thigh fat concentrates on the outer and inner portions of the contour.




Clinical Background – Thighs

■ Customized cooling cups for the thighs

- The front, inner and outer thigh regions comprise mostly of muscle, where tension is strong.
- The average cooling cup's deep suction capabilities is likely to cause pain and a bruising following procedure.
- Due to the thick muscle layer, cooling energy may not be evenly distributed and concentrated on the intended areas.
- Considering the above features of the intended target area, suitable and well-designed cooling cup is needed.


is



■ A thigh treatment that only treats the superficial layers, but takes curvatures of the contour into account.



Combinatory Treatments - Thighs



Inner Thighs

Outer Thighs & Below the Buttocks



BASELINE




4 Weeks Post Prodedure

4. Cooling Cup features

Comparison

- Thighs: Wing vs. Arc-T
 - Torso: Wing vs. Wing Plus
 - Double Chin: Wing Mini vs. Flat Mini
-

4. Thighs: Wing VS. Arc-T



Due to the anatomical features of the thighs, where the muscles are thick and broad, it is difficult to suction this area properly using Wing cup.



Post Cup Removal



Post Gel Pad Removal



Bruising may occur due to inappropriate suction using Wing Cup on thighs areas



4. Thighs: Wing VS. Arc-T

Case study: 15 subjects received Wing or Arc-T treatment and its comparison of treatment efficacy and improvement

• Wing




BASELINE



4 Weeks Post Procedure

• Arc-T



BASELINE



4 Weeks Post Procedure

Clinical Cases Comparison-



VS



Wing

Arc-T

174.5 cm ²	Treatment Area	378 cm ²	2 times ↑
90 min.	Time	45 min.	½ ↓
-0.91cm	Results	-1.15 cm	27% ↑
Bruising +++++	Side Effects	Bruising -	Near Zero
++++	Pain	+	75% ↓
4	Patient Satisfaction (1~10)	7	1.8times ↑
Improved	Patient Evaluation	Much Improved	25% ↑

4. Thighs: Wing VS. Wing plus

- Difference?



VS



Extended outer edges that can suction in more fat.



Cup Outline



Blunt-end outer edges which is suitable for smaller fat areas.

Large Fat Pockets → Abdomen

Treatment Area

Small Fat Pockets → Back

Suitable on areas on curving contours with partial fat.



Feature



Applicable for smaller areas where the vertebrae and ribs obstruct.

Cases where gel pad is vacuumed with an ineffective suction of target fat.


Recommended to use Wing Plus when cooling cup is larger than target fat.

***Note**




Appropriate for use on subjects with an overall smaller body size.

4. Double Chin: Wing mini VS. Flat mini

- Small, curving region located at the center of the submentum.



- Double chin (side) area with less curvature and small fat pockets.



5. Postures for Effective Double Chin Treatments

- Double Chin Treatment

- Safe & Comfortable Treatment Postures with Wing Mini & Flat Mini

- Wing mini Double Chin (Front)



45° Incline

- Flat mini Double Chin (Side)



30° Incline



5. Various Patient Postures During Treatment

5. Various Patient Postures During Treatment

- Wing



5. Various Patient Postures During Treatment

- Wing Plus



5. Various Patient Postures During Treatment

- Flat



5. Various Patient Postures During Treatment

- Flat mini



5. Various Patient Postures During Treatment

- Wing mini



5. Various Patient Postures During Treatment

- Arc-T



5. Various Patient Postures During Treatment

- Wing Max



Upper Abdomen



Lower Abdomen