

QUICK REFERENCE GUIDE

ATHENA



The parameters suggested here are not an absolute guide in clinical practice. The operator should be based on their own clinical experience and professional judgment to performing all processing proposed in this document.

LEN.COM.BR

ATHENA MODE				
HANDPIECE	MODE	FLUENCY	RATE OF RETREATING (mm/SHOT)	COMMENTS
ATHENA 360°	COLLIMATED	1,5 – 2,5 J/cm ²	4 - 5	Circular beam; 3.5 mm
ATHENA 90°	FRACTIONAL	30 – 40 mJ/mtz	8 - 9	SPOT 8 mm
ATHENA 90°	COLLIMATED	3 – 4 J/cm ²	8 - 9	SPOT 8 mm

ATHENA MODE				
HANDPIECE	MODE	FLUENCY	RATE OF RETREATING (mm/SHOT)	COMMENTS
ATHENA 90°	FRACTIONAL	30 – 40 mJ/mtz	8 - 9	SPOT 8 mm
ATHENA 90°	COLLIMATED	3 – 4 J/cm ²	8 - 9	SPOT 8 mm

LOCAL WHITENING AND SKIN TIGHTENING			
MODE	ENERGY	NUMBER OF SHOTS	PHOTOTYPE
FRACTIONAL INLIFT	30-40 mJ/mtz	30-80 per side	I-VI

For whitening outer lips and groin:

ENERGY	PULSE TIME	SPOT SIZE	APPLICATIONS	SHOT FREQUENCY
2-3 mJ/mtz	300 µs	8/400 mtz	1	1-3hz



ETHEREA-MX[®]
DUALMODE[®]
ATHENA[®]

VERSION 1.0 - FEBRUARY, 2017



Uma empresa
Adavium Medical