

Main Indication	Treatment Area	Intensity(%)	RF Duration(ms)	Depth(mm)	Delay Time(ms)	Pass
Scar / Rhytids	Forehead	20 - 40	100 - 120	1 - 1.5	300	2
	Nose	20 - 40	100 - 120	1 - 1.5	300	2
	Periorbital	20 - 40	100 - 120	1 - 1.5	300	2
	Cheek	40 - 60	150	1.5	300	2
	Lip	40 - 60	150	1.5	300	2
	Chin	40 - 60	150	1.5	300	2
	Neck	20 - 40	150	1.5	300	2
	Body	40 - 60	150	1.5	300	2

※ MOST IMPORTANT THING TO AVOID COMPLICATIONS: FOLLOW THE CAUTION BELOW AND DON'T MOVE THE TIP TOO FAST!!

- √ The tip head must be firmly pressed on the skin before pressing the foot pedal and until the needles are completely retracted.
- √ Use single mode and press on the foot pedal per each shoot.
- √ Never be faster than 1sec and be very careful with interval mode.
- √ Delay time must be more than twice from the value of RF duration.
- √ Depth must be adjusted to insulation of the needles: more than 1mm.
- √ If it smells burn, stop treatment or decrease applied parameters.
- √ Never use more than 60% intensity except hyperhidrosis.

※ Intensive treatment for cheek, lip, chin & body you can use; 40-60% intensity – 250ms RF duration – 2.5mm depth – 500ms delay time