



# Hollywood Package

These are suggested combination treatments from experienced clinicians. Please note that some treatment steps could be skipped according to the patient's skin response/condition.

## Hollywood Toning Package

<b>1.</b>  <b>PTP+ Premium Package</b>	<div> <div>1<sup>st</sup> session</div> <div>2<sup>nd</sup> session</div> <div>3<sup>rd</sup> session</div> </div> <p>This package improves pigment in various sizes. ▲<b>QS mode</b>: one pass to the entire face, ▲<b>Q-PTP mode</b>: second pass to the area without erythema reaction, ▲<b>Q4 mode</b>: third pass to the area without erythema reaction, ▲<b>Spectra mode</b>: lightly perform multi passes to the entire face.</p>
<b>2.</b>  <b>Melasma Package</b> [1-2 weeks]	<p>▲ Selectively destroy the melasma pigments in the dermis using <b>QS mode</b> in the early sessions, ▲ Increase clinical effect by combining <b>Q-PTP mode</b> to the remaining pigments, ▲ Switch to <b>Q3 and Q4 mode</b> to prevent melasma recurrence, ▲ <b>Spectra mode</b> to improve the overall dermal condition.</p> <p>**** <b>Note</b>: In the case of dark or melasma with a high risk of rebound, the use of QS mode could worsen the condition. In such a case, it is recommended to start the package with Q4 mode.</p>
<b>3.</b>  <b>Prevention of Recurrence Package</b> [1-2 weeks]	<p>A treatment package for melasma patients whose conditions have been improved. ▲<b>Q3 mode</b>: targets remained pigments and deliver gentle thermal energy ▲<b>Spectra mode</b>: improves the dermal condition in between the package ▲<b>MDF treatment</b>: rejuvenation treatment to improve overall skin texture.</p>
<b>4.</b>  <b>Rejuvenation Package</b> [1-2 weeks]	<p>A rejuvenation package for patients who finished melasma treatment. ▲<b>Q4 mode</b>: gentle laser toning to improve skin condition, ▲<b>MDF treatment</b>: induces collagen remodeling, improve wrinkles and enlarged pores, and uneven skin texture.</p>

## Mild Enlarged Pores treatment

<b>1.</b>  <b>MDF Rejuvenation Package</b> [1-2 weeks]	<div> <div>1<sup>st</sup> session</div> <div>2<sup>nd</sup> session</div> <div>3<sup>rd</sup> session</div> </div> <p>This is a very mild rejuvenation treatment that does not require use of anesthetic cream.</p> <ul style="list-style-type: none"> <li>• <b>Endpoint</b>: Mild Erythema</li> <li>• <b>Parameter</b>—1. (Full face): MDF hp, 1064 nm, Q4 mode, Step III, 6 mm, 3.4-4.5 J, 10 Hz, 2-3 passes</li> <li>2. (Partial area): MDF hp, 1064 nm, Q4 mode, step I, 6 mm, 3.0-3.6 J, 10 Hz, 2-3 passes</li> <li>3. (Full face): Collimated hp, 1064 nm, Spectra mode, 8 mm, 2.8 J, 10 Hz, 2-3 passes</li> </ul>
<b>2.</b>  <b>DOE Rejuvenation Package</b> [1-2 weeks]	<p>This is a very mild DOE rejuvenation treatment that does not require use of anesthetic cream.</p> <ul style="list-style-type: none"> <li>• <b>Endpoint</b>: Mild Erythema</li> <li>• <b>Parameter</b>-(Full face): DOE hp, 1064 nm, Q3 or Q4 mode, 1.0-1.2 J, 10 Hz, 2-3 passes</li> </ul>