



Case Report:

Traumatic scar treatment using ETHEREA-MX[®] with the ProDeep[®] handpiece.

Bruna Lenz, dermatologist
Belo Horizonte, MG, Brazil.

White paper case report described by dermatologist Bruna Lenz, from Belo Horizonte, Minas Gerais (MG), Brazil. The given case is of a patient complaining of traumatic scarring. Therapy comprised the use of the ETHEREA-MX[®] platform with the ProDeep[®] handpiece, with no associations before, during or after treatment. The non-ablative ProDeep[®] fractional LASER treats subdermal collagen stimulation without the inconvenient post-treatment effects of non-ablative fractionated LASERS.

Therapy:

- ▶ **1st session:** Carried out on October 26, 2017, with spot 8/100 mtz/cm² and 115 mJ/mtz with pulse time of 5 ms;
- ▶ **2nd session:** 30 days after the first, with spot 8/100 mtz/cm² and 125 mJ/mtz with time pulse of 5 ms;
- ▶ **3rd session:** Also 30 days after the second, with spot 8/100 mtz/cm² and 135 mJ/mtz with time pulse of 5 ms;
- ▶ **4th session:** A little more than 30 days after the third treatment session, with spot 8/100 mtz/cm² and 140 mJ/mtz with pulse time of 5 ms.



Pre-treatment

After the 1st treatment session.



Pre-treatment

After the 4th treatment session.



After the 2nd treatment session.

After the 3rd treatment session.