

PIGMENTARY LESION						
PHOTOTYPE	SPOT	LESION	FLUËNCY	$\lambda$	TIME OF SHOT	COOLING
I-II	40x12 mm	<i>superficial</i>	14-16 J/cm <sup>2</sup>	515 nm	10 ms	V
I-II	40x12 mm	<i>superficial</i>	13-17 J/cm <sup>2</sup>	515 nm	15 ms	V
III	40x12 mm	<i>superficial</i>	13-16 J/cm <sup>2</sup>	515 nm	15 ms	V
I-II	12X12 mm	<i>superficial</i>	14-17 J/cm <sup>2</sup>	515 nm	10 ms	V
I-II	12x12 mm	<i>superficial</i>	15-19 J/cm <sup>2</sup>	515 nm	15 ms	V
III	12x12 mm	<i>superficial</i>	14-16 J/cm <sup>2</sup>	515 nm	15 ms	V
I-II	8 mm	<i>superficial</i>	14-18 J/cm <sup>2</sup>	515 nm	10 ms	V
I-II	8 mm	<i>superficial</i>	15-19 J/cm <sup>2</sup>	515 nm	15 ms	V
III	8 mm	<i>superficial</i>	14-17 J/cm <sup>2</sup>	515 nm	15 ms	V

VASCULAR LESIONS						
PHOTOTY PE	SPOT	LESION	FLUËNCY	$\lambda$	TIME OF SHOT	COOLING
I-II	40x12 mm	< 0.3 mm	13-16 J/cm <sup>2</sup>	515 nm	10 ms	V
I-II	40x12 mm	> 0.3 mm	14-17 J/cm <sup>2</sup>	515 nm	15 ms	V
I-II	12x12 mm	< 0.3 mm	14-18 J/cm <sup>2</sup>	515 nm	10 ms	V
I-II	12x12 mm	> 0.3 mm	16-18 J/cm <sup>2</sup>	515 nm	15 ms	V
I-II	8 mm	< 0,3 mm	17-19 J/cm <sup>2</sup>	515 nm	10 ms	V
I-II	8 mm	> 0.3 mm	18-21 J/cm <sup>2</sup>	515 nm	15 ms	V
III	40x12 mm	< 0.3 mm	13-15 J/cm <sup>2</sup>	515 nm	10 ms	V
III	40x12 mm	> 0.3 mm	13-16 J/cm <sup>2</sup>	515 nm	15 ms	V
III	12x12 mm	< 0.3 mm	14-17 J/cm <sup>2</sup>	515 nm	10 ms	V
III	12x12 mm	> 0.3 mm	14-18 J/cm <sup>2</sup>	515 nm	15 ms	V
III	8 mm	< 0.3 mm	15-18 J/cm <sup>2</sup>	515 nm	10 ms	V
III	8 mm	> 0.3 mm	15-19 J/cm <sup>2</sup>	515 nm	15 ms	V