

QUICK REFERENCE GUIDELINE



These treatment guidelines are based on physician feedback and clinical experience gained in clinical studies and practical use. The information is provided as a guide only and is not prescriptive for any patient, indication or treatment. The guidelines are not designed to be a substitute for clinical training, and the system should only be operated by qualified practitioners who have received appropriate training and have thoroughly reviewed the Operator Manual that shipped with the system. When using the system, always observe Light/Energy-tissue interaction and clinical endpoints to determine appropriate settings.

LEN.COM.BR

MILD TO MODERATE WRINKLES

SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
8/100 mtz/cm ²	70-90 mJ/mtz	3-5 ms	1-2	moderate
8/100 mtz/cm ²	90-120 mJ/mtz	3-5 ms	1	moderate
8/100 mtz/cm ²	120-140 mJ/mtz	3-10 ms	1	severe

*If cooling is not used, limit the power to 100 mJ/mtz and the frequency to 1.5 Hz;

PHOTODAMAGE AND PIGMENTATION

SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
10/400 mtz/cm ²	10-16 mJ/mtz	3-5 ms	2-3	mild
10/400 mtz/cm ²	14-20 mJ/mtz	3-5 ms	2	mild
10/400 mtz/cm ²	10-14 mJ/mtz	3-10 ms	1-3	moderate

*If cooling is not used, limit the power to 15 mJ/mtz and the frequency to 1.5 Hz;

SCARS

SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
10/400 mtz/cm ²	10-16 mJ/mtz	3 ms	2	moderate
8/100 mtz/cm ²	120-150 mJ/mtz	3-5 ms	1	severe

*if cooling is not used, limit the power to 100 mJ/mtz and the frequency to 1.5 Hz; **more recent atrophic scars have better therapeutic indication;

STRETCH MARKS

SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
8/100 mtz/cm ²	90-120 mJ/mtz	3 ms	2	moderate
8/100 mtz/cm ²	120-150 mJ/mtz	3-5 ms	1	severe

*if cooling is not used, limit the power to 100 mJ/mtz and the frequency to 1.5 Hz; **more recent red stretch marks have better therapeutic indication;

MELASMA

SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
10/400 mtz/cm ²	10-16 mJ/mtz	10 ms	2	mild
10/400 mtz/cm ²	20 mJ/mtz	5-10 ms	2	mild

*whitening formulations may be combined during immediate post-treatment to improve the results;

ALOPECIA

SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
10/400 mtz/cm ²	10-16 mJ/mtz	10 ms	2	mild

*drug infiltration or topical substances may be combined (drug delivery) during immediate post-treatment to improve the results;



ETHEREA-MX[®]
PRODEEP[®]



An Advium
Medical
company