

QUICK REFERENCE GUIDELINE



These treatment guidelines are based on physician feedback and clinical experience gained in clinical studies and practical use. The information is provided as a guide only and is not prescriptive for any patient, indication or treatment. The guidelines are not designed to be a substitute for clinical training, and the system should only be operated by qualified practitioners who have received appropriate training and have thoroughly reviewed the Operator Manual that shipped with the system. When using the system, always observe Light/Energy-tissue interaction and clinical endpoints to determine appropriate settings.

LEN/COMBER

MILD TO MODERATE WRINKLES				
SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
8/100 mtz/cm ²	30-40 mJ/mtz	10-15 ms	1-2	moderate
8/100 mtz/cm ²	40-50 mJ/mtz	15 ms	1	moderate
8/100 mtz/cm ²	50-60 mJ/mtz	15 ms	1	severe

PHOTODAMAGE AND PIGMENTATION				
SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
10/400 mtz/cm ²	10-15 mJ/mtz	10-15 ms	2-3	mild
10/400 mtz/cm ²	4-12 mJ/mtz	10-15 ms	2	mild
10/400 mtz/cm ²	6-14 mJ/mtz	10-15 ms	1-3	moderate

SCARS				
SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
8/100 mtz/cm ²	40-50 mJ/mtz	15 ms	1	moderate
8/100 mtz/cm ²	50-60 mJ/mtz	15 ms	1	severe

STRETCH MARKS				
SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
8/100 mtz/cm ²	40-50 mJ/mtz	15 ms	1	moderate
8/100 mtz/cm ²	50-60 mJ/mtz	15 ms	1	severe

MELASMA				
SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
10/400 mtz/cm ²	3-5 mJ/mtz	15 ms	2	mild
10/400 mtz/cm ²	4-6 mJ/mtz	15 ms	1-2	mild

ALOPECIA				
SPOT SIZE	POWER	PULSE WIDTH	PASSES	EDEMA AND ERYTHEMA
10/400 mtz/cm ²	3-5 mJ/mtz	15 ms	2	mild



ETHEREA-MX[®]
GO-SMOOTH[®]

VERSION 1.2 - NOVEMBER 2017



An Adavium
Medical
company