


VARIOUS NEEDLE CARTRIDGES



MTR-AC-25

25pin non-insulated
micro-needle cartridges

Treatment Area ; 10*10(mm)

Max Output ; 25watt

Main Usage ;
mostly main indication



MTR-AC-C25

25pin semi-insulated
micro-needle cartridges

Treatment Area ; 10*10(mm)

Max Output ; 25watt

Main Usage ;
Required for higher thermal
effect at dermis/subcutaneous
tissues




MTR-AC-10

10pin non-insulated
micro-needle cartridges

Treatment Area ; 0.25*18(mm)

Max Output ; 25watt

Main Usage ;
stretch mark, naso labial folds,
crow's feet



MTR-AC-C64

64pin semi-insulated
micro-needle cartridges

Treatment Area ; 17.5*17.5(mm)

Max Output ; 70watt

Main Usage ;
Required for higher thermal
effect at dermis/subcutaneous
tissues