

FREQUENTLY ASKED QUESTIONS

What is Nano-Ultrasound?

Sonocare emits tiny vibrations of ultrasound waves that provide a massaging effect within the cell, this causes the cells to vibrate against each other. The thermal energy created as a result is delivered deep within the skin allowing for a no downtime, comfortable, soothing treatment with amazing results.

What is Multi-Wave Technology?

Sonocare's unique Multi-Wave Technology optimises the benefits of Nano-Ultrasound by delivering the energy deeper into the tissue for greater effect. Sonocare's treatments alternate between two frequencies for deep tissue tightening in the jowl area as well as superficial skin tightening for around the eyes.

Treatment time and interval of Sonocare treatment?

Each treatment takes approximately 30 minutes. 1-2 treatments per week are usually recommended. Your clinician will discuss targeted treatments options for enhanced results.

Any downtime?

Sonocare is a breakthrough in skin hydration technology. Treatment feels like a hot stone massage. You are able to get back to normal life right after the treatment.

When can I see the results?

You can see the results immediately after treatment – specifically for lymphatic drainage/swelling reduction, tightening and hydration. For some applications, it may need more sessions. Normally, 8-10 times treatment is recommended and continuous treatment maximise treatment effect.

ULTRASOUND TECHNOLOGY DESIGNED TO RESTORE SKIN QUALITY

- Tightening without dehydration
- Moisturising without injection
- Superb results without side-effects
- Short treatment without downtime
- Broken capillaries & bruising
- Wrinkle reduction

Clinic Details

Distributed by:

Cryomed
Aesthetics

cryomed.com.au

Please consult your medical practitioner to see if you are suitable for treatment. Potential risks and side effects need to be discussed with your medical practitioner before undergoing treatment with Sonocare.

© Cryomed Aesthetics Australia 2020


sonocare
NANO ULTRASOUND

SKIN TIGHTENING AND MOISTURISING



**TIGHTEN AND MOISTURISE
YOUR SKIN
SEE THE BENEFITS OF
NEW NANO ULTRASOUND**

SEE & FEEL THE DIFFERENCE

STEP 1

After the Sonocare hydration gel is applied to the skin, treatment is commenced. You will experience a warm sensation of ultrasound upon treatment commencement. Sonocare is a synergistic treatment utilising two probes for combined treatment effects.

STEP 2

The first ultrasound probe provides a non thermal effect which is glided along the skin surface stimulating lymphatic flow to reduce swelling and offers a soothing 'recovery effect' for the skin. During application of the second probe you will experience mild heating for fibroblast stimulation through this thermal effect.

STEP 3

Fibroblasts are important for skin regeneration for dry, rough, and dehydrated skin. Immediate effect on the skin is noticed in fine lines, tightening, flushing and hydration. Long term effects include collagen remodelling which improves skin elasticity for sagging of double chin, fine lines and overall skin tightening.

In these circumstances (among others), Sonocare treatment might need to be postponed. You should seek medical advice if any of the following apply to you:

- Pregnancy or lactation
- Open skin or wound
- Active skin infection or rash
- Pacemaker
- Some vascular disease causing excessive bleeding
- Some blood disorders or medications

CUTTING EDGE TECHNOLOGY

Sonocare is the world's leading dual frequency ultrasound technology designed to restore skin quality, and also delivers treatment benefits for stubborn fat, pigmentation and redness.


Sonocare utilises 10MHz of nano-ultrasound to deliver thermal energy to the epidermis and dermal layer with Mutli Wave Technology. Thermal energy at 45~48°C is delivered by nano-vibration where skin renewal occurs without damaging the skin surface. Sonocare's unique Multi-Wave Technology emits two frequencies alternatively from one transducer for highly convenient treatments.

ADVANTAGE OF 10MHZ

- 10MHz allows more targeted treatment in the epidermis and upper dermis layers for hydration 3MHz for deep dermis targeted treatment for tightening and soothing.
- 10MHz+3MHz Combination for lymphatic drainage requirements for optimised results.
- More cavitation caused in smaller cells.

Ideal depth for collagen generation

- Create heat at 0.1~0.5mm depth.
- Safer since the heat comes from inside the cell.




FOR REAL PATIENTS



Before and after skin tightening



Before and after skin tightening



Before and after soothing



Before and after body tightening