

Multiple treatments in just one machine. Combination of radio frequency and ultrasound.





Combination of treatments

the combination of radio frequency and ultrasound technologies in the dalyance equipment puts vanguard and innovation at the service of its clients with an unbeatable profitability.

Proved results

the key to a perfect body is firmness, a well-defined contour and a soft, luminous skin. dalyance treatments guarantee these results producing a pleasant sensation of confidence in its patients.

Quality guaranteed

dalyance has been designed and manufactured by a team of highly qualified professionals. The criteria used for the development of this equipment is based on the experience of the **cocoon medical** brand in the development of technology and noninvasive medical equipment.

Multiple treatments in just one machine. Combination of radio frequency and ultrasound.



Firms up, smooths and defines.

Visible results, our best guarantee

Two technologies in just one machine

dalyance equipment combines radio frequency and ultrasound in just one machine.

The purpose of radio frequency is to improve the quality of the skin and stimulate collagenesis by heating the adipose tissue and tissue at cutaneous level through heat diffusion.

Cavitation is a treatment which uses low frequency ultrasound. Micro bubbles are formed accumulating a great quantity of energy. These are able to break the membrane of the adipose cell on making contact with the adipocytes, without harming the adjacent structures and respecting the integrity of the blood and lymphatic systems, transforming the fat from a solid to a liquid state to be eliminated by the lymphatic system through the urinary tract.

How does dalyance work?

dalyance is an integral cosmetic solution which combines radio frequency and cavitation technologies for body contouring and facial rejuvenation treatments.

The combination of both technologies allows heating of the area to be treated by radio frequency while the ultrasound applies tridimensional pressure on the adipose tissue, obtaining a much more effective effect on the treated area.

Facial rejuvenation: dalyance treatment is the most effective for a truly rejuvenated skin. It performs an anti-aging process allowing the patient to notice results from the first session. The effects on the skin are visible in:

- Firmer and rejuvenated skin
- Facial contouring
- Firming up of the skin
- Reduction of the facial contour
- Firming up of the face and wrinkle reduction

Body contouring: dalyance treatments stand out for being noninvasive and painless. Its innovating technology allows unbeatable results to be obtained, making it possible to observe the loss from the first session. Furthermore dalyance allows the treatment of cellulite in all its stages: hard, flaccid and oedematous, obtaining the best result for each patient.

Body treatments offered by dalyance are:

- Cellulite reduction
- Firming up
- Improved aspect of the skin
- Body contouring

Intuitive interface

cocoon medical, in their commitment to vanguard and design, have developed very intuitive and easy-to-use interfaces, helping users to obtain the best results in their treatments.

dalyance has a very intuitive interface which facilitates the work of the user, reducing training time and helping to obtain maximum performance of the machine.



Very intuitive and easy-to-use interfaces for their equipment which help the user to obtain the best results from their treatments.

Tendency, innovation and design

Multiple accessories for an integral treatment

dalyance has been designed to offer both facial and body treatments. Thus, it works with different accesories and interchangeable heads, making work easier and adapting itself to the individual needs of each client.

Accesories for facial rejuvenation: the equipment has a facial hand piece with a interchangeable 35mm head (monopolar and bipolar). The resuls offered are:

- Redefinition of face contour
- Reduction of wrinkles and expression lines
- Reduction of bags and rings
- Lifting of perioocular region
- Filling of wrinkles in nasolabial region

Body 6 pin multipolar: allows the destruction of the adipose tissue and helps eliminate localized fat. Its use is recommended in the middle treatment sessions.



Monopolar: this head helps to firm up the skin obtaining a smooth and elastic effect. Its use is recommended in the final sessions.



Monopolaire

Accessories for body contouring: dalyance offers body contouring in 3 phases: drainage, reduction and firming up, allowing the number of sessions to be adapted to the client's expectations.

Body 4 pin multipolar: suitable for draining the tissue and favouring the elimination of liquids. Its use is recommended in the first sessions.





Concave ultrasound handpiece: it works by destroying the adipocyte and favouring the reduction of adipose tissue. It is indicated to treat areas with higuer tissue.



Flat ultrasound hand piece: it works by destroying the adipocyte and favouring the reduction of adipose tissue. It is indicated to treat areas with less tissue.



Results

Unbeatable results

dalyance has been designed by a team of highly qualified professionals working to make the achievement of results the company motto.

cocoon medical has a unit for training and clinical monitoring where the results of the level of therapeutic response are checked in vivo, together with the application of protocols to obtain the best results.

The response obtained from dalyance treatments is a gradual improvement in the treated area, visible from the first session and long lasting once finished. The combination of the dalyance technologies favour wrinkle reduction and body contouring.

dalyance works for facial rejuvenation, highlighting the following: reduction of wrinkles in the periocular region, filling of nasogenian wrinkles, reduction of fine lines such as barcode, eyelid lift and tightening of areas such as neck and neckline.

As for body contouring, dalyance can work in three phases: drainage, reduction and firming up. It obtains an integral recontouring effect and an unbeatable result.

Before

After



43 years old woman. Soft cellulite and sagging. After 4 sessions.



58 years old woman. Soft cellulite and sagging. After 4 sessions.

Before







43 years old woman. Hard cellulite. After 4 sessions.

Before

Before



After



72 years old man. Deep wrinkles and cutaneous elastosis grade IV. After 5 sessions.

Technical specifications

Basic operation

Type of equipment:

ultrasound and radio frequency equipment

Technical characteristics

Maximum power Radio frequency power Ultrasound power Ultrasound frequency Radio frequency frequency 300 W 180 W 120W 36 KHz - 40 KHz 1MHz

Power supply

power supply Marking 00 -240 Vac 50 - 60 Hz CE

Directives

Directive 2006/95/CEE (LVD) on Low Voltage Directive 2004/108/CEE (EMC) on Electromagnetic Compatibility Directive 2002/95/CEE (RoHS) Directive 2002/96/CEE (WEEE)

Physical characteristics

Materials Equipment weight Dimensions Screen metal and methacrylate finishing 20 kg. 430 x 400 x 300 mm touch screen 10,4"



make your patients' dreams come true

cocoon medical is a Spanish company, based in Barcelona, created with the aim of offering technology and innovation to the service of cosmetic medicine, capable of manufacturing non-invasive medical equipment to satisfy the needs of the most demanding clients.

Built-up by a team of highly qualified professionals, **cocoon medical** has become a leading worldwide company due to its continued ability for innovation and development. We place our maximum trust in our partners, both national and international, in order to get to understand the needs of each market, offering the most innovating solutions to our customers.

cocoon medical has the highest national and international quality certifications, offering the best service in their field.

www.cocoonmedical.com



High Technology Products PASSATGE MASOLIVER, 24-28 08005 BARCELONA, SPAIN

Customer Service: T. +34 934 588 566 info@cocoonmedical.com



Follow us on:

f Y in 🛅