

SECRET SmartCure MICRONEEDLING RF PARAMETERS

NEEDLE	INDICATION	NEEDLE SIZE	NEUTRAL ELECTRODE REQUIRED	INTENSITY	RF	DISTANCE B/N	TOPICAL REQD	REPEAT	# SESSIONS
K3i	Vascular <150 microns	0.075mm	NO	20	30ms	3mm	NO	3 - 4 weeks	2 - 3
AC-1	Acne	1.5mm	YES	20-30	Continuous		YES	3 - 4 weeks	1 - 2
	Periorbital Wrinkles	1.5mm	YES	20	Continuous	1mm	YES	4 weeks	2 - 3
	Perioral Wrinkles	1.5mm	YES	10	Continuous	1mm	YES	4 - 8 weeks	1 - 3
Subcision	Scars	19G, 22G, 27G	YES	10	Continuous	Fan	YES	4 weeks	2 - 3