



Hollywood Package

Cryomed
aesthetics

These are suggested combination treatments from experienced clinicians. Please note that some treatment steps could be skipped according to the patient's skin response/condition.

Hollywood Toning Package

1. PTP+ Premium Package	<div> <div>1st session</div> <div>2nd session</div> <div>3rd session</div> </div> <p>This package improves pigment in various sizes. ▲ QS mode: one pass to the entire face, ▲ Q-PTP mode: second pass to the area without erythema reaction, ▲ Q4 mode: third pass to the area without erythema reaction, ▲ Spectra mode: lightly perform multi passes to the entire face.</p>
2. Melasma Package [1-2 weeks]	<p>▲ Selectively destroy the melasma pigments in the dermis using QS mode in the early sessions, ▲ Increase clinical effect by combining Q-PTP mode to the remaining pigments, ▲ Switch to Q3 and Q4 mode to prevent melasma recurrence, ▲ Spectra mode to improve the overall dermal condition.</p> <p>**** Note: In the case of dark or melasma with a high risk of rebound, the use of QS mode could worsen the condition. In such a case, it is recommended to start the package with Q4 mode.</p>
3. Prevention of Recurrence Package [1-2 weeks]	<p>A treatment package for melasma patients whose conditions have been improved. ▲ Q3 mode: targets remained pigments and deliver gentle thermal energy ▲ Spectra mode: improves the dermal condition in between the package ▲ MDF treatment: rejuvenation treatment to improve overall skin texture.</p>
4. Rejuvenation Package [1-2 weeks]	<p>A rejuvenation package for patients who finished melasma treatment. ▲ Q4 mode: gentle laser toning to improve skin condition, ▲ MDF treatment: induces collagen remodeling, improve wrinkles and enlarged pores, and uneven skin texture.</p>

Mild Enlarged Pores treatment

1. MDF Rejuvenation Package [1-2 weeks]	<div> <div>1st session</div> <div>2nd session</div> <div>3rd session</div> </div> <p>This is a very mild rejuvenation treatment that does not require use of anesthetic cream.</p> <ul style="list-style-type: none"> • Endpoint: Mild Erythema • Parameter—1. (Full face): MDF hp, 1064 nm, Q4 mode, Step III, 6 mm, 3.4-4.5 J, 10 Hz, 2-3 passes 2. (Partial area): MDF hp, 1064 nm, Q4 mode, step I, 6 mm, 3.0-3.6 J, 10 Hz, 2-3 passes 3. (Full face): Collimated hp, 1064 nm, Spectra mode, 8 mm, 2.8 J, 10 Hz, 2-3 passes
2. DOE Rejuvenation Package [1-2 weeks]	<p>This is a very mild DOE rejuvenation treatment that does not require use of anesthetic cream.</p> <ul style="list-style-type: none"> • Endpoint: Mild Erythema • Parameter-(Full face): DOE hp, 1064 nm, Q3 or Q4 mode, 1.0-1.2 J, 10 Hz, 2-3 passes