

# CLATUU $\alpha$

## Double Chin treatment Guide



# Patients selection

- Make sure it's not just sagging skin but adipose tissue. It is most suitable for the procedure of fat tissue over 2cm thick.
- Check the health information such as injection or skin disease, heart disease, and cooling allergy before the procedure. Identify the site of treatment desired by the patient.

# Select the right cup for patient

Wing Mini



Flat Mini



# Marking the area



- Each grid is used to select the most suitable cup. Generally, the flat mini cup is suitable for both sides of the jaw, and the Wing mini cup is suitable for the center of the lower jaw.
- The appropriate cup is selected and designed in consideration of the shape of the patient's double jaw or bending.
- Draw the outer part of the grid with a dark pen to mark the cup position.

# Set the parameters

## Wing Mini



Suction	Cooling	Time
2~3	4~5	40 Minutes

## Flat Mini



Suction	Cooling	Time
2~3	4~5	40 Minutes



Suction	Cooling	Time
2~3	4~5	40 Minutes

# Apply the gel pad(CLATUU shield 360 C type)



- Check gel pad C type and open it.
- Spread the gel pad on the treatment area. When attaching, make sure that the design part is located at the center.
- Remove the mesh between the gel pads and the air layer between the skin and the gel pads.

# Place the Cushion – Wing mini



# Place the Cushion – Flat mini





# Apply the cup on the chin



# Placement of the Mini cup

Lateral Jaw Procedure



Lateral Jaw Procedure



# Massage the area



- After removing the gel pad, use a warm towel to wipe the gel off your skin. Massage the area gently for 2 to 3 minutes by hand.

# Massage the area



# Common immediate side effects

- redness
- firmness
- pinching sensation
- minor swelling
- tenderness
- bruising
- stinging
- tingling
- numbness
- skin sensitivity
- mild to moderate pain
- muscle cramping

# Additional rare complications that may occur after the treatment.

- pain that starts several days after the treatment and usually resolves when numbness goes away, or after about a month
- scarring or changes to skin color
- reduced movement or response of the tongue following procedures involving the neck or chin
- lower lip muscle weakness following procedures involving the neck or chin
- reduced saliva production and dry mouth following procedures involving the chin or neck