

CuRAS MLA HP Treatment Guide

This guide is a reference data (mainly SKIN TYPE III) so need to be adjusted according to the skin condition and response

	Indication	Mode	Spot size(mm)	Fluence (J/ cm ²)	Frequency (Hz)	Treatment Guide
Acne Scar & Large pore	Toning	Q1064	7mm	1.15 ~1.23 J/cm ²	10 Hz	<ul style="list-style-type: none"> ✓ Applied to relatively thick skin with pigmentation caused by trouble acne ✓ Approximately 2,000 to 2,500 shots at full face
	Skin Texture	Q1064	7mm	1.00 J/ cm ²	15 Hz	<ul style="list-style-type: none"> ✓ Partial treatment in areas where Erythema has not occurred. (Endpoint ; Erythema) ✓ Deeply delivers low level energy repetitively by 4 ~ 5 stacking on deep scars and pores
	Acne Scar, Pore	Q1064	7mm	2.51~2.59 J/cm ²	5 Hz	<ul style="list-style-type: none"> ✓ Delivers the energy repeatedly by 2 ~ 3 stacking (or 2pass) to the depressed scars or pores. The endpoint is when the Petechia is occurring immediately after ; a small (1 – 2 mm) red or purple spot on the skin ✓ In the case of deep scarred particularly, additional 2 passes treatment is required till check the Edema immediately after as the endpoint
Melasma Blemish & Pigment	Toning	Q1064	7mm	0.84 ~ 0.92 J/cm ²	10 Hz	<ul style="list-style-type: none"> ✓ When treating a spot where a freckle has occurred, perform the procedure by lowering of the energy level than other areas. For example, if the default value is 0.92 J/cm², the freckle's treatment energy is proper with 0.84 J/cm² ✓ Approximately 2,000 to 2,500 shots at full face
	Melasma	PTP1064	7mm	0.7 J/cm ²	10 Hz	<ul style="list-style-type: none"> ✓ Selectively 4~5 overlap to the areas where freckle, blemishes, and pigmentations are occurring.
	Melasma, Blemish	PTP1064	7mm	1.00 J/cm ²	10 Hz	<ul style="list-style-type: none"> ✓ If the spot is dark than other area, apply 1 ~ 2 pass lightly around the certain spot area. ✓ The endpoint is Erythema ✓ For Synergistic effect, add 1 ~ 2pass lightly rejuvenation (M1064) at the finishing step.
Rejuvenation & Wrinkle	Toning	Q1064	7mm	0.84 ~ 0.92 J/cm ²	10 Hz	<ul style="list-style-type: none"> ✓ Combine Q1064 (approx. 2,000 shots) for toning effect and M1064 (approx .1,000 shots) for skin regeneration, elasticity and wrinkle improvement
	Rejuvenation	M1064	7mm	1.19 ~1.79 J/cm ²	10	<ul style="list-style-type: none"> ✓ After the pigment treatment, lightly add 1 ~ 2 passes at the finishing step for synergy effect. ✓ The endpoint is Erythema
	<p>*Care is advised for dermal pigment including Melasma. Proper assessment and diagnosis of Melasma is important for best clinical outcome. Consider Zoom hand piece for moderate to severe Melasma.</p>					