



Indications: Professional and Amateur Tattoos. Always start tattoo removal treatments with the 1064 - treat the entire tattoo including all colors before using the other wavelengths on the tattoo.

Spectra Tattoo Removal Checklist / Log

Wavelength (nm)	1064	532	585 (Gold HP) Use 532 wavelength mode	650 (Red HP) Use 532 wavelength mode
Spot Size	**8 – 4 mm	2.6 mm	2mm	2mm
Fluence (J/cm2)	1.6 – 2.2j	.28 – .90j	4.0j	3.0j
Treatments	6 – 10	3 – 4	2 – 4	2 – 4
Color	Black/Blue, Brown & all colors	Purple, Red/Orange	Yellow, Sky Blue	Green

Name: _____ Skin Type: I II III IV Provider: _____

Consent: _____ Allergies: _____ Topical: _____ Local: _____

Photo: _____ Area Treated: _____ Tattoo Colors: _____

TX #	Date	Topical On / Off	Test Spot	Wave-length	Spot Size	Fluence	Pulse Count	Provider / Laser Notes
1			<input type="checkbox"/>					
			<input type="checkbox"/>					
2			<input type="checkbox"/>					
			<input type="checkbox"/>					
3			<input type="checkbox"/>					
			<input type="checkbox"/>					
4			<input type="checkbox"/>					
			<input type="checkbox"/>					
5			<input type="checkbox"/>					
			<input type="checkbox"/>					
6			<input type="checkbox"/>					
			<input type="checkbox"/>					
7			<input type="checkbox"/>					
			<input type="checkbox"/>					
8			<input type="checkbox"/>					
			<input type="checkbox"/>					
9			<input type="checkbox"/>					
			<input type="checkbox"/>					
10			<input type="checkbox"/>					
			<input type="checkbox"/>					