

FILMED
LABORATOIRES

AUSTRALIA



NANOSOFT™

A NEW EXPERIENCE

IN SKIN REJUVENATION

DELICATE ZONES DIFFICULT TO TREAT

THE MAIN WOMEN'S CONCERNS

YOUR SKIN REFLECTING YOU



**SOME SKIN ZONES ARE
MORE FRAGILE
THAN OTHER**

**48% DARK CIRCLES
35% CROW'S FEET**

**35% LINES AROUND
THE MOUTH AND LIPS**

**27% NECK
(FIRMNESS &
ELASTICITY)**

EYE CONTOUR - PERIORAL AREA - NECK & DECOLLETAGE



First areas
to show
ageing signs



Continuous
movement



Difficult to target
efficiently

NANOSOFT™ DEVICE

NEW EXCLUSIVE PENETRATION TECHNOLOGY

NANOSOFT™ BY FILLMED, A REAL REVOLUTION IN AESTHETICS

EASY TO HANDLE
Blue line for injection's
orientation

OPTIMAL DESIGN
Can be mounted on any
standard syringe



CONTROLLED & SECURE
INTRADERMAL DELIVERY
3 microneedles, 0.6mm in
length

Certified medical device

Dispenser of 30 sterile devices Sterilization by EtO (Ethylene-Oxide)

NANOSOFT™, THE SMALLEST NEEDLE ON THE MARKET


3 TIMES SMALLER
THAN AN INJECTION
NEEDLE



* Arcane study 2018 - 8,646 women (aged 18-75).

UNIQUE MICRO-NEEDLING TECHNOLOGY


CHRONO-DIFFUSION FROM THE PAPULAS



1/ Into the upper dermis



2/ Creation of a papula



3/ Chrono-diffusion

4/ The papulas disappear after
24h max

NANOSOFT™

BETTER RESULTS IN LESS TIME

INTRADERMAL WITH NANOSOFT™

PERIORBITAL AREA



by Dr. Valérie PHILIPPON

NECK



by Dr. Riekie SMIT



by Dr. Riekie SMIT



by Dr. Riekie SMIT

BEFORE

AFTER 1 SESSION

VISIBLE RESULTS
EVEN AFTER
1 SESSION

HOW TO USE NANOSOFT™



1/ The blue line must always face you.



2/ Stretch the skin and insert the device at a 45° angle.



3/ Inject slowly and expect strong resistance, a tiny papula will form. Repeat every centimetre to cover the full area.

TRAINING MOVIE



BEAUTY IS AN ART, WE MADE IT A SCIENCE

FILLMED

A DIVISION OF THE CRYOMED GROUP IN AUSTRALIA AND NEW ZEALAND

1300 346 448 E: orders@cryomed.com.au W: fillmed.com.au