

FREQUENTLY ASKED QUESTIONS

HOW DOES DUALMODE TREATMENT FEEL?

Laser skin rejuvenation treatments involve only very mild discomfort. You may experience a prickling sensation as the skin tissue is heated. In some cases, an anaesthetic cream may be applied prior to treatment.

WHAT HAPPENS DURING TREATMENT?

The operator will move the handpiece carefully over the treatment area. You can expect this to be repeated in a number of passes.

HOW MANY TREATMENTS WILL I NEED?

The condition of your skin will dictate how many treatments you need to undergo. Between three and five is an average number.

HOW SOON WILL I SEE RESULTS?

You will see improvement within a few weeks, but your skin will continue to renew itself and will actually look even better a month or so after your last treatment. For this reason, repeat treatments should be spaced several weeks apart.

For best ongoing results, you should always protect your skin from direct sun exposure. Always cover up and use a good SPF sunblock when out in the sun.

GIVE YOURSELF THE SKIN YOU DESERVE WITH DUALMODE

ER:YAG 2940nm LASER
FOR NO DOWNTIME SKIN REJUVENATION

- Suits all skin types
- Low downtime
- Minimal discomfort
- Continuous improvement
- Low risk of side effects

Clinic details:

Cryomed
Aesthetics

cryomed.com.au

Please consult your practitioner to see if you are suitable for treatment. Potential risks and side effects need to be discussed with your practitioner before undergoing treatment with DualMode.

© Cryomed Aesthetics 2016

REJUVENATE YOUR SKIN



**LASER SKIN
REJUVENATION**
IMPROVES TONE,
TEXTURE & EARLY
SIGNS OF AGEING



024CV1

RESTORE YOUTHFULNESS

Your skin is the front line of defence against a hostile environment. The result: sun damage, uneven tone, superficial lines and wrinkles. It's time to help your skin feel young again.

Unlock the potential for rejuvenation within your own skin cells. DualMode's technology works to stimulate your skin to act like the youthful skin you used to have.

In the past, aggressive treatments such as chemical peels and ablative procedures involved extensive downtime and discomfort. They were not suitable for all skin types.

DualMode rejuvenation treatments harness your body's own potential to heal itself. They help break down the damaged cells so your body can remove them, then they stimulate your body's natural process to produce fresh, elastic cells that support a more youthful skin.

STEP 1

Your skin cells are stimulated with a precisely controlled energy source. This triggers your skin cells to renew themselves, producing quantities of the chemicals that gave your skin its youthful elasticity and fullness.

STEP 2

Your body replaces the damaged skin cells with new, healthy cells.

STEP 3

The result is younger, firmer, clearer-looking skin. Whether you need a more even skin tone, smoothing of superficial lines and wrinkles, or firming – or all three – the DualMode can turn back the clock for your skin.

THE DUALMODE ER:YAG 2940NM LASER FOR SKIN REJUVENATION



Before



After



Before



After



Before



After

RESULTS YOU CAN SEE AND FEEL



Before



After



Before



After

** Individual results may vary based on the aggressiveness of treatment and skin condition.*